

An LED Approach for Measuring the Photocatalytic Breakdown of Dye on Titanium Dioxide Coated Surfaces

Robert E. Ryan,¹ Lauren W. Underwood,² Duane O'Neal,² Mary Pagnutti,¹ and Bruce A. Davis³
¹Innovative Imaging and Research, Bldg. 1103, John C. Stennis Space Center, MS 39525
²Science Systems and Applications, Inc., Bldg. 1105, John C. Stennis Space Center, MS 39529
³Department of Homeland Security, Science and Technology Directorate, Washington, DC 20528



ABSTRACT

A simple method to estimate the photocatalytic reactivity performance of spray-on titanium dioxide (TiO₂) coatings for transmissive glass surfaces was developed. This novel technique provides a standardized method to evaluate the efficiency of photocatalytic material systems under a variety of illumination levels. To date, photocatalysis assessments have generally been conducted using mercury black light lamps. Illumination levels for these types of lamps are difficult to vary, consequently limiting their use for assessing material performance under a diverse range of simulated environmental conditions. This new technique uses a stable ultraviolet (UV) gallium nitride (GaN) light emitting diode (LED) array instead of a traditional black light lamp to initiate and sustain photocatalytic breakdown. This method was tested with a UV-resistant dye (Crystal Violet) applied to a TiO₂ coated glass slide. Experimental control is accomplished by applying Crystal Violet to both TiO₂ coated slides and uncoated control slides. A slide is illuminated from the dye side by the UV LED array at various light levels representative of outdoor and indoor conditions. To monitor degradation of the dye over time, a temperature-stabilized white light LED, whose emission spectrum overlaps with the dye absorption spectrum, is used to illuminate the opposite side of the slide. Using a spectrometer, the amount of light from the white light LED transmitted through the slide as the dye degrades is monitored as a function of wavelength and time. In this way, the rate of dye degradation for photocatalytically coated versus uncoated slide surfaces can be compared. Results demonstrate that the dye transmission increased much more rapidly on the photocatalytically coated slides than on the control uncoated slides, and that dye degradation is dependent on illumination level. For photocatalytic activity assessment purposes, this experimental configuration and methodology minimizes many external variable effects and enables small changes in absorption to be measured. This research demonstrates the use of this alternative, innovative LED light source design over traditional mercury black light lamp systems and non-LED lamp approaches. This novel technology begins to address the growing need for a standard method that can assess the performance of photocatalytic materials before deployment for large-scale, real world use.

INTRODUCTION TO PHOTOCATALYSIS

In the presence of UV light (320–400 nm) and water, the photocatalyst TiO₂ has two unique qualities:

- Breakdown Properties**—accelerated decomposition of compounds:
 - Organic and inorganic compounds (i.e., Volatile Organic Compounds-VOCs, Polychlorinated Biphenyls-PCBs, pesticides, Nitrogen Oxides-NO_x)
 - Antimicrobial/antibactericidal (i.e., bacterium, fungi, viruses)
- Self-Cleaning Properties**—creates highly washable surface state:
 - Contact angle between the photocatalyst surface and water is reduced
 - Induces superhydrophilicity




Figure 1. Photocatalysis

Reduction-oxidation reaction occurs at the surface of the catalyst and generates:
-OH, hydroxyl radicals and
O₂⁻, superoxide ions
These molecules break down the structure of the microorganisms and of inorganic and organic (VOCs) compounds and convert them into CO₂ and H₂O.

PROBLEM AND PROPOSED SOLUTION

Problem: Limited data available on performance of commercially available photocatalytic materials under a variety of environmental conditions.

Solution:

- Time series transmission spectra of dye breakdown are taken to determine photocatalytic reactivity.
- Crystal Violet dye used to assess degradation capability (use of other dyes may also be possible).
- Traditional UV simulator black light lamp replaced with a GaN UV LED array to allow for variation in light level.

Figure 2. Transmission Through Dye

